

ELECTROSURGERY COLPOSCOPY

Colour coded
Lietz
loops

- Ultra- fine rigid tungsten wire for enhanced control
- Nonstick balls minimise adhesion



COLPOSCOPY LOOPS - *easy to identify the size by colour*



PMS 1742
Loop electrode green
10mm x 10mm loop
Pack of 25 - single use



PMS 1746
Loop electrode orange
15mm x 8mm loop
Pack of 25 - single use



PMS 1747
Loop electrode pink
15mm x 12mm loop
Pack of 25 - single use



PMS 1743
Loop electrode blue
15mm x 15mm loop
Pack of 25 - single use



PMS 1748
Loop electrode green
20mm x 10mm loop
Pack of 25 - single use




PMS 1751
Loop electrode grey
20mm x 12mm loop
Pack of 25 - single use




PMS 1752
Loop electrode yellow
20mm x 13mm loop
Pack of 25 - single use



PMS 1745
Loop electrode red
20mm x 15mm loop
Pack of 25 - single use



PMS 1744
Loop electrode beige
20mm x 20mm loop
Pack of 25 - single use



PMS 1735

Square loop electrode
Light blue
5mm x 5mm
Pack of 25 - single use

PMS 1732


Square loop electrode
Light blue
10mm x 8mm
Pack of 25 - single use



Failsafe identification



Crush resistant packaging



Easy to secure fingergrrips

nonstick BALL ELECTRODES



PMS 2908 **nonstick** ball electrode
3mm ball, 110mm long. Pack of 10 - single use



PMS 2907 **nonstick** ball electrode
4mm ball, 110mm long. Pack of 10 - single use



PMS 2906 **nonstick** ball electrode
5mm ball, 110mm long. Pack of 10 - single use



PMS 2908L **nonstick** ball electrode
3mm ball, 130mm long. Pack of 10 - single use



PMS 2907L **nonstick** ball electrode
4mm ball, 130mm long. Pack of 10 - single use



PMS 2906L **nonstick** ball electrode
5mm ball, 130mm long. Pack of 10 - single use

STANDARD BALLS




PMS 1761 ball electrode
3mm ball, 110mm long. Pack of 10 - single use



PMS 1762 ball electrode
5mm Ball, 110mm long. Pack of 10 - single use


The importance of nonstick balls



Conventional electrode

Stainless steel electrodes

Conventional electrodes have a tendency for localised heating which may result in tissue becoming attached and subsequently breaking off once the electrode is removed causing post operative bleeding.



Prima Nonstick electrode

Nonstick electrodes

- Cleaner
- Faster
- Safer

Prima
Inspiring innovation

Tel: +44 (0) 1453 752626

Fax: +44 (0) 1453 756750

Email: sales@prima-medical.com

Web: www.prima-medical.com

Prima Medical Limited,
Mill 16, Bath Road Trading Estate,
Stroud, Gloucestershire,
GL5 3QF, UK

Distributed by: