



*progress in surgery*

Titanium  
ARTERIAL & VENOUS  
BULLDOG  
CLAMPS



LAPAROSCOPY

# TITANIUM ARTERIAL & VENOUS BULLDOG CLAMPS

## Advantages:

- Different pressure color coding (A-Artery / V-Vein)
- Fast recognition and differentiation, confusion is avoided
- Lightweight
- Non-magnetic
- Highest biocompatibility
- Mechanical strength
- Corrosion-resistance
- DeBaKey serration

# A



Titanium Bulldog Clamp, straight, Code A ● pressure: 475-500gr.

|               |          |
|---------------|----------|
| 300-100-773TA | 17mm jaw |
|---------------|----------|

|               |          |
|---------------|----------|
| 300-100-774TA | 25mm jaw |
|---------------|----------|

|               |          |
|---------------|----------|
| 300-100-771TA | 45mm jaw |
|---------------|----------|

|               |          |
|---------------|----------|
| 300-100-772TA | 60mm jaw |
|---------------|----------|

Not for sale in Europe

# V



Titanium Bulldog Clamp, straight, Code V ● pressure: 300-375gr.

|               |          |
|---------------|----------|
| 300-100-773TV | 17mm jaw |
|---------------|----------|

|               |          |
|---------------|----------|
| 300-100-774TV | 25mm jaw |
|---------------|----------|

|               |          |
|---------------|----------|
| 300-100-771TV | 45mm jaw |
|---------------|----------|

|               |          |
|---------------|----------|
| 300-100-772TV | 60mm jaw |
|---------------|----------|

Not for sale in Europe



300-101-775

300-410-775

## Endoscopic Bulldog Clamp Applier

|             |   |
|-------------|---|
| 300-101-775 | ∅ 10mm, 330mm, rotatable, with Luer Lock, serrated jaws |
|-------------|---|

|             |                  |
|-------------|------------------|
| 300-410-775 | same, detachable |
|-------------|------------------|

|             |                               |
|-------------|-------------------------------|
| 300-401-775 | insert only (for 300-410-775) |
|-------------|-------------------------------|

## ENDOSCOPIC Klein™ BULLDOG CLAMPS

### Advantages:

- Secure holding of occluded tissue structure
- Occlusion of vessel and other tissue without additional port
- Atraumatic cross clamping of the colon
- DeBakey serration



Vascular Bulldog Clamp, straight

|             |          |
|-------------|----------|
| 300-100-773 | 17mm jaw |
|-------------|----------|

|             |          |
|-------------|----------|
| 300-100-774 | 25mm jaw |
|-------------|----------|



Bulldog Clamp for bowel surgery, straight

|             |          |
|-------------|----------|
| 300-100-771 | 45mm jaw |
|-------------|----------|

|             |          |
|-------------|----------|
| 300-100-772 | 60mm jaw |
|-------------|----------|



Vascular Bulldog Clamp, curved

|             |          |
|-------------|----------|
| 300-100-776 | 17mm jaw |
|-------------|----------|

|             |          |
|-------------|----------|
| 300-100-777 | 25mm jaw |
|-------------|----------|



Bulldog Clamp for bowel surgery, curved

|             |          |
|-------------|----------|
| 300-100-778 | 45mm jaw |
|-------------|----------|

|             |          |
|-------------|----------|
| 300-100-779 | 60mm jaw |
|-------------|----------|



Satinsky Bulldog Clamp, for use with 12.5mm trocar system

|             |          |
|-------------|----------|
| 300-100-786 | 60mm jaw |
|-------------|----------|



Satinsky Bulldog Clamp, for use with 12.5mm trocar system

|             |          |
|-------------|----------|
| 300-100-789 | 60mm jaw |
|-------------|----------|



*progress in surgery*

## RZ Medizintechnik GmbH

Unter Hasslen 20  
D-78532 Tuttlingen

Phone: +49 (0)7462/9470-0  
Fax: +49 (0)7462/9470-50

[sales@rz-medizintechnik.com](mailto:sales@rz-medizintechnik.com)  
[www.rz-medizintechnik.com](http://www.rz-medizintechnik.com)

Media Code Number: 021610